



# QualyLive

Live with more quality

# QualyBone

Dental



# 100%

TRICALCIUM  
PHOSPHATE (BETA-TCP)

■ POROUS SYNTHETIC CERAMIC DESIGNED FOR BONE CAVITY FILLING

■ INDUCES REGENERATION AND BONE GROWTH, STIMULATING THE PROLIFERATION AND DIFFERENTIATION OF OSTEOBLASTS

■ DUE TO ITS MESH CONFIGURATION IT HAS A HIGH POROSITY FACILITATING THE PROLIFERATION OF THE BONE CELLS, INTO THE VOIDS

■ THIS FORMAT GIVES IT EXCELLENT MECHANICAL RESISTANCE

■ CELL ADHESION OBSERVED 5 DAYS AFTER BONE GRAFTING

## INDICATIONS

- FILLING DEFECTS OR BONE CAVITIES
- VOLUME INCREASE IN BONE DEFECTS THAT ARE NOT INTRINSIC TO THE STABILITY OF THE BONE STRUCTURE
  - Reconstruction of tumour voids and cyst defects
  - Crestal augmentation
  - Alveolar regeneration
  - Regeneration of periodontal defects
  - Sinus lift





# WHY CHOOSE QUALYBONE TCP?

## SAFETY

100% SYNTHETIC AND RESORBABLE

## RESORBABLE

FROM 6 TO 12 MONTHS, DURING THIS PERIOD IT IS REPLACED BY NEW BONE TISSUE

## RADIOPAQUE

EXCELLENT RADIOLOGICAL CONTROL FEATURES, DURING THE BONE REGENERATION PROCESS AND UNTIL ITS TOTAL OSTEOINTEGRATION

## HIGH COHESIVITY AND VOLUME CONSERVATION

BEING VERY HYDROPHILIC GIVES IT AN EXCELLENT PARTICLE COHESIVENESS, WHICH IN TURN ALLOWS FOR THE INITIAL CAVITY VOLUME CONSERVATION


## VASCULARIZATION

INDUCES AN EFFICIENT VASCULARIZATION PROCESS SOON AFTER THE INTERVENTION

## DIRECT APLICATION

NO MEMBRANE REQUIRED. DUE TO ITS PHYSICAL AND MECHANICAL PROPERTIES

## EASY TO USE



GEOMETRY	SIZE	QUANTITIES	REF.
GRANULES	0,1 - 0,5 mm	0,5g	QBT010505G
		0,5 x 5	QBT010505P
		1g	QBT010510G
		1g x 5	QBT010510P
GRANULES	0,5 - 1 mm	0,5g	QBT050105G
		0,5g x 5	QBT050105P
		1g	QBT050110G
		1g x 5	QBT050110P
BLOCKS	5x10x15 mm	0,5 g	QBT010205G
		0,5g x 5	QBT010205P
		1g	QBT010210G
		1g x 5	QBT010210P
		1 UNIT	QBT051015B



Qualy Live  
 Rua Soares dos Reis, 1183 - R/C Eq  
 4430-241 - Mafamude - Vila Nova de Gaia  
 Porto - Portugal  
 ☎ +351 215861750  
 ☎ +351 935652996  
 ✉ [qualylive@qualylive.com](mailto:qualylive@qualylive.com)  
 🌐 [www.qualylive.com](http://www.qualylive.com)