



QualityLive

Live with more quality

QualityBone

Dental



75%

HYDROXYAPATITE

25%

TRICALCIUM
PHOSPHATE (BETA-TCP)

■ POROUS SYNTHETIC CERAMIC DESIGNED FOR BONE CAVITY FILLING

■ INDUCES REGENERATION AND BONE GROWTH, STIMULATING THE PROLIFERATION AND DIFFERENTIATION OF OSTEOBLASTS

■ DUE TO ITS MESH CONFIGURATION IT HAS A HIGH POROSITY FACILITATING THE PROLIFERATION OF THE BONE CELLS, INTO THE VOIDS

■ THIS FORMAT GIVES IT EXCELLENT MECHANICAL RESISTANCE

■ CELL ADHESION OBSERVED 5 DAYS AFTER BONE GRAFTING

INDICATIONS

- FILLING DEFECTS OR BONE CAVITIES
- VOLUME INCREASE IN BONE DEFECTS THAT ARE NOT INTRINSIC TO THE STABILITY OF THE BONE STRUCTURE
 - Reconstruction of tumour voids and cyst defects
 - Crestal augmentation
 - Alveolar regeneration
 - Regeneration of periodontal defects
 - Sinus lift



WHY CHOOSE QUALYBONE BCP?

SAFETY

100% SYNTHETIC AND RESORBABLE

RESORBABLE

FROM 12 TO 24 MONTHS, DURING THIS PERIOD IT IS REPLACED BY NEW BONE TISSUE

RADIOPAQUE

EXCELLENT RADIOLOGICAL CONTROL FEATURES, DURING THE BONE REGENERATION PROCESS AND UNTIL ITS TOTAL OSTEOINTEGRATION

HIGH COHESIVITY AND VOLUME CONSERVATION

BEING VERY HYDROPHILIC GIVES IT AN EXCELLENT PARTICLE COHESIVENESS, WHICH IN TURN ALLOWS FOR THE INITIAL CAVITY VOLUME CONSERVATION


VASCULARIZATION

INDUCES AN EFFICIENT VASCULARIZATION PROCESS SOON AFTER THE INTERVENTION

DIRECT APLICATION

NO MEMBRANE REQUIRED. DUE TO ITS PHYSICAL AND MECHANICAL PROPERTIES

EASY TO USE



GEOMETRY	SIZE	QUANTITIES	REF.
GRANULES	0,1 - 0,5 mm	0,5g	QBB010505G
		0,5 x 5	QBB010505P
		1g	QBB010510G
		1g x 5	QBB010510P
GRANULES	0,5 - 1 mm	0,5g	QBB050105G
		0,5g x 5	QBB050105P
		1g	QBB050110G
		1g x 5	QBB050110P
BLOCKS	5x10x15 mm	0,5 g	QBB010205G
		0,5g x 5	QBB010205P
		1g	QBB010210G
		1g x 5	QBB010210P
		1 UNIT	QBB051015B



Qualy Live
 Rua Soares dos Reis, 1183 - R/C Eq
 4430-241 - Mafamude - Vila Nova de Gaia
 Porto - Portugal
 ☎ +351 215861750
 ☎ +351 935652996
 ✉ qualylive@qualylive.com
 🌐 www.qualylive.com

Ver. EN_Dental_BCP_09062022